

## **Aniara Diagnostica Launches [hit-assay.Com](http://hit-assay.com) Providing Specific Products for Detection of Heparin Induced Thrombocytopenia**

*Aniara Diagnostica is pleased to announce a comprehensive line of standardized ELISA Methods optimized for measuring Heparin- dependent antibodies of the IgG and IgGAM isotypes in human plasma.*

WEST CHESTER, Pa. ([PRWEB](http://PRWEB)) November 23, 2018 -- Aniara Diagnostica, a multi-national distributor of diagnostic and research agents, as well as screening tools, has strengthened its commitment to meeting today's highest standards for product quality and cost-efficient delivery. Recently the company announced a new site [hit-assay.com](http://hit-assay.com) which highlights the availability of standardized ELISA Methods for the detection of Heparin-dependent antibodies.

### **ZYMUTEST™ HIA IgG**

Enzyme Immuno-assay (HIT) designed for measuring Heparin-dependent antibodies of the IgG isotype, in human plasma.

- ZYMUTEST™ HIA IgG (96 tests) (ARK040A)
- ZYMUTEST™ HIA MonoStrip IgG (32 tests) (ARK041A)

### **ZYMUTEST™ HIA IgGAM**

The ZYMUTEST™ HIA IgGAM kit is a screening assay that measures globally all human heparin-dependent antibodies of the IgG, IgM or IgA isotypes.

- ZYMUTEST™ HIA IgGAM (96 tests) (ARK040D)
- ZYMUTEST™ HIA MonoStrip IgGAM (32 tests) (ARK041D)

Aniara distributes Hyphen BioMed's entire line of diagnostic and research reagents for Thrombosis, Hemostasis, and Autoimmunity.

The entire Aniara product line is available for purchase at <http://www.aniara.com>. The Company distributes products in North America, Central America, South America, Sweden, Norway, Denmark, Finland, Lithuania, Latvia, Estonia and Iceland.



**Contact Information**

**Maria Waldron**

Aniara Diagnostica

<http://www.aniara.com>

9054840680

**Online Web 2.0 Version**

You can read the online version of this press release [here](#).