

Product Spotlight: Aniara Diagnostica would like to announce the availability of Protein L Ligand Leakage ELISA kit

Aniara Diagnostica is pleased to announce the availability of Protein L Ligand Leakage ELISA kit for the detection and quantification of Protein L in Immunoglobulin (Ig) or Ig-fragment containing solutions.

WEST CHESTER, Ohio ([PRWEB](#)) February 24, 2021 -- Aniara Diagnostica, a multi-national distributor of diagnostic and research agents, as well as screening tools, has strengthened its commitment to meeting today's highest standards for product quality and cost-efficient delivery.

For Research Use Only. Not For Use In Diagnostic Procedures.

Aniara Diagnostica is now the exclusive distributor of Protein L kits in the USA and Canada.

The Protein L Ligand Leakage ELISA Kit has been designed to detect and quantify Protein L in Immunoglobulin (Ig) or Ig-fragment containing solutions. It is a sandwich ELISA with microtiter strips coated with an affinity purified anti-Protein L IgY-antibody. The Protein L Ligand Leakage ELISA kit comes in two versions – the standard kit (A10-0027) and the XL kit (A10-0028-1).

A10-0027 Protein L Ligand Leakage ELISA (Standard Kit) 12 strips (96 wells) [click here for more information](#)

A10-0028-1 Protein L Ligand Leakage ELISA XL (1 kit) 12 strips (96 wells). The XL kit contains five vials of Protein L reference instead of one and may be used for process adaptation purposes. [click here for more information](#)

The entire Aniara product line is available for purchase at <http://www.aniara.com>. The Company distributes products in North America, Central America, South America, Sweden, Norway, Denmark, Finland, Lithuania, Latvia, Estonia and Iceland.



Contact Information

Maria Waldron

Aniara Diagnostica

<http://www.aniara.com>

5133423018

Online Web 2.0 Version

You can read the online version of this press release [here](#).