

## **Product Spotlight: Aniara Diagnostica would like to announce the availability of Recombinant Hirudin**

*Aniara Diagnostica is pleased to announce the availability of Recombinant Hirudin, a thrombin inhibitor.*

WEST CHESTER, Ohio ([PRWEB](#)) November 12, 2019 -- Aniara Diagnostica, a multi-national distributor of diagnostic and research agents, as well as screening tools, has strengthened its commitment to meeting today's highest standards for product quality and cost-efficient delivery.

For Research Use Only. Not for Use in Diagnostic Procedures.

Recombinant Hirudin is a Thrombin inhibitor, which can be used in any in vitro laboratory research application where an efficient inhibition of Thrombin is required.

ARE120A Recombinant Hirudin (1 x 16,000 ATU) [click here for more information](#)

ARE120B Recombinant Hirudin (1 x 160,000 ATU) [click here for more information](#)

The entire Aniara product line is available for purchase at <http://www.aniara.com>. The Company distributes products in North America, Central America, South America, Sweden, Norway, Denmark, Finland, Lithuania, Latvia, Estonia and Iceland.



**Contact Information**

**Maria Waldron**

Aniara Diagnostica

<http://www.aniara.com>

5133423018

**Online Web 2.0 Version**

You can read the online version of this press release [here](#).