

## **Product Spotlight: Aniara Diagnostica would like to announce the availability of T cell ELISPOT assays**

*Aniara Diagnostica is pleased to announce the availability of T cell ELISPOT assays for COVID-19 vaccine research.*

WEST CHESTER, Ohio ([PRWEB](#)) June 11, 2020 -- Aniara Diagnostica, a multi-national distributor of diagnostic and research agents, as well as screening tools, has strengthened its commitment to meeting today's highest standards for product quality and cost-efficient delivery.

For Research Use Only. Not for Use in Diagnostic Procedures.

T cell ELISPOT assays are extremely effective for COVID-19 vaccine research. Because Human and monkey research is critical in the COVID-19 field, IFN-gamma T cell ELISPOT assays are an attractive tool to determine the number of COVID-specific IFN-gamma producing T cells in the circulation of human and non-human primates.

A cytokine assay that has also been reported to be important for COVID-19 research is the IL-6 ELISA (cytokine ELISA kit).

ACT230-PR5 Human IFN- $\gamma$  ELISPOT kit [click here for more information](#)

ACT121-PR5 Monkey IFN- $\gamma$  ELISPOT kit [click here for more information](#)

ACT640-20 Human IFN- IFN- $\gamma$  ELISPOT antibody pair [click here for more information](#)

ACT605-20 Monkey IFN- $\gamma$  ELISPOT Antibody Pair [click here for more information](#)

ACT205A Human IL-6 ELISA [click here for more information](#)

The entire Aniara product line is available for purchase at <http://www.aniara.com>. The Company distributes products in North America, Central America, South America, Sweden, Norway, Denmark, Finland, Lithuania, Latvia, Estonia and Iceland.



**Contact Information**

**Maria Waldron**

Aniara Diagnostica

<http://www.aniara.com>

513-342-3018

**Online Web 2.0 Version**

You can read the online version of this press release [here](#).